

# SE-1010 PC ECG



Wired

Wireless

## SE-1010 PC ECG



### About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Edan Instruments, Inc. | 3/F-B, Nanshan Medical Equipments Park,  
1019# Nanhai Rd., Shenzhen | 518067 P.R. China  
+86.755.26898326 | www.edan.com.cn | info@edan.com.cn



ENG-ECG-PC ECG-V3.0-20150120

© 2015 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission.






# SE-1010

## PC ECG

With a data acquisition module and a set of comprehensive software program for acquiring, analyzing, storing and reviewing 12-channel ECGs, SE-1010 PC ECG system provides all-in-one solutions for both resting ECG and stress test.



Wired and wireless solutions make your options more flexible.



Multi formats ECG reports export and can be sent by email to streamline your workflow.



### Ability to work with both PC & Tablet

EDAN's PC ECG offers you flexibility to view ECG traces on both PC and Tablet. **(PC: Windows 7 32/64-bit, Windows Vista 32/64-bit, Windows XP SP2/SP3. Tablet: Windows 7 32/64-bit)**

### Wired Solution

Wired solution provides you with simplified installation.

### Wireless Solution

Wireless solution provides the freedom of patients and clinicians during the stress test.

### 12-channel Resting ECG

- Clear and reliable ECG waveforms acquisition
- Accurate measurements and interpretation tested with CSE/MIT/AHA database

### Comprehensive Analysis

- Vector ECG, Timed VCG, Signal Average Electrocardiogram (SAECG)
- Heart Rate Turbulence (HRT), QT interval Dispersion
- High Frequency ECG, Frequency Spectrum cardiogram


### Versatile Report

- PDF, FDA-XML, SCP, DICOM, JPG, BMP and DOC format reports transferrable via E-mail
- Compatible with GDT systems

### Stress test (Optional)

- Wired and wireless solutions
- Compatible with reputable treadmills and bicycles
- Support exercise blood pressure monitor
- Customized exercise protocol in addition to Bruce and Modified Bruce

PC-based ECG SE-1010 offers you an affordable and portable 12-channel ECG solution with complete data analysis functions.



QT Interval Dispersion

HF ECG

FCG

HRV

SAECG