

SonoTrax Vascular

Ultrasonic Pocket Doppler



EDAN

A world of potential

SonoTrax Vascular Ultrasonic Pocket Doppler

SonoTrax Vascular has been widely adopted as standard equipment for the assessment of peripheral vascular disease. The versatile, compact and portable SonoTrax Vascular allows fast and simple but effective assessment of both arterial and venous disease.

Features

- Easy-to-use
- Compact design
- Support 1.5 V AA battery and Rechargeable Ni-MH Battery
- Accurate blood flow detection with high-fidelity sound
- Support earphone output
- Support 4 MHz, 5 MHz, 8 MHz waterproof probes
- Interchangeable waterproof probe (IPX8)

Specifications

- Audio Output Power: 1 W
- Size: 34 (D) x 89 (W) x 141 (H) mm
- Weight: about 300 g (including the battery)

