

MFM-CMS

Central Monitoring System



MFM-CMS
Central Monitoring System



About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



A world of potential

Edan Instruments, Inc. | 3/F-B, Nanshan Medical Equipments Park,
1019# Nanhai Rd., Shenzhen | 518067 P.R. China
+86.755.26898326 | www.edan.com.cn | info@edan.com.cn

© 2015 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission.



ENG-PM-MFM-CMS-V2.0-20150323



A world of potential



- Dual-display monitoring
- User-friendly interface and optimized workflow
- Connection via LAN/WLAN to 128 bedside monitors
- Support HL7/XML
- Monitoring data statistic and analysis
- Drug dose, Hemodynamic, Oxygenation, Renal function, Ventilation calculation
- Bi-directional communication with EDAN bedside patient monitors
- Time synchronization
- Remote NIBP measurement
- Batch management
- Patient info setup
- Alarm setup

MFM-CMS

Central Monitoring System

MFM-CMS brings you a comprehensive central monitoring solution for all EDAN bedside monitors. It offers seamless monitoring management from transportation to hospitalization, from intensive care to general care, bringing easy access to well-organized monitoring technology.

Storage capacity for each patient:

720 NIBP Review	720 Alarm Event	240 h Trend	72 h Full-disclosure Waveform	12 h Dynamic Short Trend
---------------------------	---------------------------	-----------------------	---	------------------------------------

Calculation Function

Five kinds of calculations are introduced in MFM-CMS to provide an overall clinical guidance including Drug Dose, Hemodynamic, Oxygenation, Renal Function and Ventilation calculation.



Statistic and Analysis

Besides monitoring statistic, the analysis for alarm event, ECG ARR and historical data can be showed in management window either in graphic or in tabular view.

